


Maintaining Your CodaAir® Systems




Changing Your Filters Regularly Is Extremely Important

4-Stage Filter System (4 Canisters) – Change Filters Every 6 Months



Semi-Annual Kit
Change 6 months after annual kit.
(4 Canisters, 2" & 4" prefilters)



Annual Kit
Change 6 months after semi-annual kit.
(4 Canisters, 2" & 4" prefilters, 6" HEPA)

Positive Pressure Unit



Catalogue #: 115V 230V	CAPP-115 CAPP-230
Semi-Annual Kit	CPPS-001
Annual Kit	CPPA-001

ECO™ 1500



Catalogue #: 115V 230V	GECO-115 GECO-230
Semi-Annual Kit	GESA-1500
Annual Kit	GEAK-1500

ECO™ 1200



Catalogue #: 115V 230V	GELT-115 GELT-230
Semi-Annual Kit	GESA-1200
Annual Kit	GEAK-1200


CodaAir® 900



The CodaAir® 900 offers you two choices of filter systems. – Choice B

Catalogue #: 24V Universal Adapter	A900-009
Semi-Annual Kit	CCSK-900
Annual Kit	CCFK-900

Shirt Box Filter System – Change Filters Every 6 Months



Semi-Annual Kit – Sold in 2 packs
Change every 6 months.
(2" pleated filter, 2" carbon/potassium filter, 4" HEPA)

CodaAir® 900



The CodaAir® 900 offers you two choices of filter systems. – Choice A

Catalogue #: 24V Universal Adapter	A900-009
Initial Filter	Included
Semi-Annual Kit	CAFK-900 (2pk filter kit)

CodaAir® 800



Catalogue #: 24V Universal Adapter	A800-008
Initial Filter	Included
Semi-Annual Kit	CAFK-800 (2pk filter kit)

Aero® 700



Catalogue #: 24V Universal Adapter	A700-000
Initial Filter	Included
Semi-Annual Kit	GCAF-702 (2pk filter kit)

Round Filters – Change Filters Every 3 Months



Quarterly Kit

Change every 3 mths.
(Sold in 1 or 2 pks)

Aero® 500



Catalogue #: 24V Universal Adapter	GCAU-000
Initial Filter	Included
Quarterly Kit – 1 pk	GCAF-001 (1pk filter kit)
Quarterly Kit – 2 pk	GCAF-002 (2pk filter kit)

Coda® Filters – Change Filters Every 30 days

Coda®2 Unit

Place on the top shelf of the incubator.



Catalogue #: Unit Only	C2IU-010
Catalogue #: Kit 1 Coda®2 Unit, 6 Coda® Filters, 1 Coda® Xtra Inline® Filter, 1 Connector Set, 1 Coda® Inline® Filter Magnet, 1 Coda® Passthru Device	C2KT-107

Coda® Filters Change every 30 days.	
Single Filter	GCIF-001
6 pk of Filters	GCIF-006
2 pk of Filters	GCIF-012



Coda® Unit

Place on the top shelf of the incubator.

Catalogue #:	GCIU-010
---------------------	----------



Coda® SP

Mount into the access hole of the incubator.

Catalogue #:	GCSP-010
---------------------	----------



Coda® Inline® Filters – Change Coda Xtra® Inline® Filters Every 3 Months – Change Coda Regular Inline® Filters Every 2 Months

Coda® Xtra Inline® Filters – Green with AldaSorb®

Attach to your incoming gas lines.

For best results, change filters every 3 months or when gas cylinders are replaced (can be used up to 6 months).

Single Filter	CXGR-001
6 pk of Filters	CXGR-006
2 pk of Filters	CXGR-012



Coda® Xtra Inline® Filters – Blue

Attach to your incoming gas lines.

For best results, change filters every 3 months or when gas cylinders are replaced (can be used up to 6 months).

Single Filter	CXBL-001
6 pk of Filters	CXBL-006
2 pk of Filters	CXBL-012



Coda® Regular Inline® Filters – Purple

Attach to your incoming gas lines.

For best results, change filters every 2 months or when gas cylinders are replaced (can be used up to 3 months).

Single Filter	CIPR-001
6 pk of Filters	CIPR-006
2 pk of Filters	CIPR-012

