



Rocket® Craft™ DUOVAC illustrated

- **HIGHLY RELIABLE:** Rocket Craft™ Pumps have a proven 30 year history of reliability and excellent value.
- **SIMPLY, SAFE OPERATION** with air-operated foot switches on a solid moulded base to replace the traditional chromed style. Lighter weight connection set is now standard across the range to improve ease of handling and reliability.
- **LARGER, CLEARER GAUGE** with increased resolution improves setting accuracy and long term reliability.
- **ULTRA QUIET** diaphragm pump makes all Rocket pumps easy to live with in the OR and typically only needs bi-annual servicing.
- **FULLY VARIABLE VACUUM** from 0-400mmHg<sup>-1</sup>. High Vacuum: 440mmHg<sup>-1</sup>
- **SINGLE USE WATER TRAP or FILTER SETS** provide the simplest and most cost effective solution to the need to reprocess reusable water traps and silicone tubing.



**LabIVF®**  
Touching Lives  
WORLDWIDE

Description		Order Code
<b>Rocket® Craft™ Suction Pump</b> complete with air operated vacuum footswitch, 5 single patient use water trap sets, mains supply cord, User Instruction Manual	220-240VAC 50Hz	<b>R29655</b>
	110-120VAC 60Hz	<b>R29654</b>
<b>Rocket® Craft™ DUOVAC Suction Pump</b> complete with air operated high/low vacuum footswitch, 5 single patient use water trap sets, mains supply cord, User Instruction Manual	220-240VAC	<b>R29660</b>
	110-120VAC	<b>R29661</b>
<b>Rocket® Single Use Water Trap Set</b> 2.5m patient connection tube with male luer connector. Supplied sterile, for single use in cartons 10 units.		<b>R57685</b>